

# SAN DIEGO MESA COLLEGE AUTOCLAVE SPORE TESTING RESULTS

## To Be Performed Monthly!

Manufacturer: \_\_\_\_\_ Model # \_\_\_\_\_

District Equipment #: \_\_\_\_\_ Location \_\_\_\_\_

Date of Test: \_\_\_\_\_ Length of Test Cycle: \_\_\_\_\_

Materials Used in Test Cycle: \_\_\_\_\_

Maximum Temperature of Autoclave Cycle: \_\_\_\_\_

### **RESULTS OF INCUBATION:**

**Follow manufacturer instructions on proper testing method and procedures:**

Temperature of incubator: \_\_\_\_\_

	Sterilized Ampule	Control Ampule
24 hours		
48 Hours		
72 hours		

**Final Results** \_\_\_\_\_ (Record as Pass or Fail)

### Key to Results:

- + indicates growth (turbidity) and yellow color. Sterilization was not accomplished. When recording state color net to (+ or -) symbol
- indicates no growth and the ample remained purple (or red) in color. When recording state color next to the (+ or -) symbol.

**Retain record for a minimum of 3 years.**